

GE MAC 1600 EKG Machine – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The GE MAC 1600 is a sophisticated ECG system equipped with the Marquette 12SL analysis program for accurate diagnosis and informed clinical decision-making. The device supports resting ECG, rhythm monitoring, arrhythmia detection, and optional exercise testing capabilities, allowing clinicians to access and review ECG data at any time.

The MAC 1600 features a full-color WVGA display showing all 12 ECG leads simultaneously, with dedicated lead keys for easy navigation between 6- and 12-lead review. Its compact design and efficient operation make it suitable for small practices and mid-sized clinics.

The system includes a new patient cable with banana lead wires, ECG clips, new battery, new power cord, thermal paper, and user manual. An optional rolling cart is available for enhanced mobility and convenience.

FEATURES

- Slide-out paper drawer can be easily loaded and quickly adjusted for US letter and A4 paper sizes.
- Advanced Lithium Ion battery technology permits three hours of continuous use.
- Dedicated lead keys permit scrolling between 6- and 12-lead review.
- Full-color WVGA display of all 12 ECG leads.
- F1-6 keys directly correspond to detail in boxes above.
- Integrated handle makes transportation easier.

SPECIFICATIONS

- Type: 7" TFT LCD module, widescreen
- Display Resolution: 800 x 480 pixels
- Height: 15 in/400 mm
- Width: 13.8 in/350 mm
- Depth: 6 in/150 mm
- Weight: Approximately 6.0 kg, with battery installed
- Voltage: 100 to 240 VAC, +10, -15%
- 12 lead simultaneous analysis

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)