

GE MAC 2000 EKG Machine – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Refurbished GE MAC 2000 is powered by the proven 12-lead Marquette Analysis Program, delivering fast and accurate ECG interpretation. Furthermore, it provides high-quality resting ECG diagnostics and offers optional stress (strain) ECG capabilities.

The MAC 2000 is designed to boost clinical efficiency and support accurate, timely assessments across various healthcare settings. The GE MAC 2000 is ideal for private practices, busy clinics, and pharmaceutical clinical trials.

Additionally, A user-friendly interface and full keyboard simplify operation for staff. Furthermore, its easy setup reduces training time, helping physicians and nurses use it confidently and efficiently.

FEATURES

- Marquette 12SL ECG Analysis Program
- PDF and XML Export Capabilities
- Large Easy Read 7 Inch, Color Display
- Up-To 200 ECG's Internal Storage
- Rechargeable Battery for Portability

SPECIFICATIONS

- Weight: 11 Lbs
- Product Dimensions: 13 Inches Width x 15.5 Inches Length x 8 Inches Height
- Paper Size: 8.5 x 11, 8.25 x 11.6, 8.43 x 11 Inches
- Paper Type: Z-fold Thermal Pre-Printed Grid, Perforation
- Printer Speed: 5, 12.5, 25, 50 Millimeters Per Second

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)