

UMF 4010 Head-Centric Procedure Chair



Stock Photo - Options May Vary

FACTORY NEW

The UMF 4010 is a head-centric procedure chair designed for clinical precision and patient comfort across a wide range of medical procedures. It features fully powered hi-lo, backrest, tilt, and leg extension functions that enable effortless positioning and improved ergonomic access for providers. The intuitive hand control with OneTouch Patient Positioning allows caregivers to adjust to preset or custom positions.

The chair includes a double-articulating headrest that ensures precise head and neck support, adapting to different procedures and patient needs. A removable, seamless upholstered top features Prefixx Protective Finish, offering durability, stain resistance, and easy cleaning to reduce patient turnover time and support infection control protocols.

Designed for functionality and durability, the UMF 4010 delivers comfort, control, and reliability for modern clinical environments while protecting patients from infection, facilities from damage, and staff from injury.

FEATURES

- Chair designed for head-centric practices
- Double-articulating headrest with height adjustability and interchangeable headrest options
- 4 power functions: hi-lo, back, tilt, leg lift & extension
- Hand control with OneTouch, patient positioning system, 2 programmable positions and auto-return
- Power leg lift for patient comfort
- Trendelenberg, Lithotomy & Semi-Fowler positions
- Removable, seamless upholstery with Prefixx protective finish
- Hospital-grade electrical receptacle with child safety cover (patient left)

SPECIFICATIONS

- Product Type: Procedure Chair
- Height: Min. 18.5" (± .5") at Seat Max. 37"
- Seat Length: 28"
- Back Length: 24.5"
- Length: 70" - 80" with Foot Section Extended
- Leg Rest Dimensions: 28" (L) x 13" (W)
- Patient Weight Capacity: 650 lbs
- Electrical: 120 VAC, 60HZ, 5 AMP

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)