

UMF Proglide4010 Procedure Chair



Stock Photo - Options May Vary

FACTORY NEW

The UMF Proglide4010 Procedure Chair is designed to create a comfortable, accessible, and hygienic environment for patients and providers. Its compact footprint, slim back cushion, and low 18.5" height meet ADA guidelines and allow easy wheelchair access.

The chair supports up to 650 lbs and features a durable steel base for enhanced stability and long-term reliability. A seamless, removable upholstered cushion supports infection control and allows for quick cleaning between patients.

The Proglide4010 Features: Dual electro-mechanical actuators power Hi-Lo, Back, Tilt, Leg Lift, and Extension adjustments with smooth, efficient movement. Two programmable positions and a return-to-home function increase workflow efficiency and reduce setup time.

FEATURES

- OneTouch WheelBase System with central-locking caster system allows chair to be moved around room for cleaning
- Double-articulating headrest supports patients for head-centric procedures
- Headrest raises and lowers and is interchangeable with a round, u-shaped, and rectangular option (U-shaped and rectangular headrest not included)
- Powered hi-lo, back, tilt, and leg extension are all controlled with a hand remote or a foot remote
- OneTouch Patient Positioning® System allows 2 positions to be programmed into memory and an auto-return

SPECIFICATIONS

- Height: Min. 18.5" (+ .5") 46.9cm (+1.27cm) at seat Max. 37" at seat (93.9cm)
- Seat Width: 28" (71.12cm)
- Back Width: 24.5" (62.23cm)
- Length: 70" - 80" with foot section extended (177.8cm — 203.2cm)
- Leg Rest: 28" W x 13" D (71.12cm W x 33.02cm D)

**Brochure and images are a general overview of the product. Actual product options may vary.*