

UMF FusionFREEDOM™ Exam Chair



Stock Photo - Options May Vary

FACTORY NEW

The UMF FusionFREEDOM™ Exam Chair meets ADA standards recommended by the U.S. Access Board for inclusive clinical care. A low seat height of 17 inches enables safer, easier, and more independent transfers for patients with mobility challenges, accommodating 95% of wheelchair users capable of self-transfer.

Many exam rooms still use non-adjustable, fixed-height tables despite over 20 million adults facing mobility limitations. FusionFREEDOM provides a practical, ADA-compliant alternative for modern healthcare settings. TruComfort™ Support Rails provide added stability, reducing fall risk and supporting patient and caregiver safety during exams.

This cost-effective solution promotes equitable, dignified care for every patient while protecting patients from infection, facilities from damage, and staff from injury.

FEATURES

- Power hi-lo
- 450lb (204kg) patient weight capacity
- 17" (43.2cm) low-access ADA seat height
- TruComfort™ Support Rails
- SmoothGlide Leg™ Extension with Privacy Pull™ Loops

SPECIFICATIONS

- Height: min 17" (43.18cm)
- Width: 28" (71.1 cm)
- Length: 58"-76" (147.3 cm-193 cm) with leg rest extended
- Weight capacity: 500 lbs. (226.8kg)
- Electrical: 120 VAC 10.5 AMP / 220 VAC 5.8 AMP, 50-60Hz

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)