

# UMF 4070 Power100 Power Exam Table



Stock Photo - Options May Vary

## FACTORY NEW

The UMF 4070 Power100 Power Exam Table combines durability, ADA accessibility, and smooth power operation for efficient patient care. Its dual electro-mechanical actuators power Hi-Lo and back adjustments for quiet, consistent positioning. A low 18.5" seat height meets ADA standards and allows safer, easier patient transfers.

Furthermore, with a 650 lb weight capacity, it supports a wide range of patients with confidence. A 7° pelvic tilt improves positioning for exams and procedures. UMF 4070 Power100 Features include a built-in drain pan, foot control, and a front drawer for convenient storage.

Removable, seamless upholstery ensures easy cleaning and better infection control. Adjustable stirrups and a front leg cushion add flexibility and comfort.

## FEATURES

- Power Hi-Lo with Air Spring Assisted Manual Backrest
- 650lb (294.84kg) Patient Weight Capacity
- Low Access ADA Height 18.5" ( ±0.5") 46.9cm (±1.27cm)
- One-function Foot Control - Power Hi-Lo
- Removable Seamless Upholstered Cushion

## SPECIFICATIONS

- Height: Min. 18.5" (± .5") 46.9cm (±1.27cm) at Seat Max. 37" at seat (93.9cm)
- Width: 28" (71.12cm)
- Length: 57" - 76" with foot section extended (144.7cm - 193cm)
- Leg Rest: 19" W x 13" D (48.26cm W x 33cm D)
- Front Storage Drawer: 19" W x 12" D (48.26cm W x 30.48cm D)

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)