

UMF 4070 Power200 Power Exam Table



Stock Photo - Options May Vary

FACTORY NEW

The UMF 4070 Power200 Power Exam Table reflects clinical research and provider input. It enhances patient care and operational efficiency. In addition, it meets ADA and U.S.

Access Board standards, supporting accessibility for a wide range of patient needs. The seamless, removable upholstered top with PreFixx® finish allows deep cleaning and easy replacement for infection prevention. Furthermore, Powder-coated steel surfaces clean easily and resist wear, promoting hygiene in high-traffic clinical settings.

The Power200 Exam Table enhances efficiency with OneTouch Patient Positioning®, offering two programmable positions and return-to-home functionality. Power Hi-Lo and a manual backrest enable smooth adjustments, while 7° pelvic tilt and stainless steel drain pan enhance precision.

FEATURES

- Power Hi-Lo with Air Spring Assisted Manual Backrest
- 650lb (294.84kg) Patient Weight Capacity
- Low Access ADA Height 18.5" (±0.5") 46.9cm (±1.27cm)
- Hand Control with OneTouch Patient Positioning® 2 Programmable Positions & Return-to-Home
- One-function Foot Control - Power Hi-Lo

SPECIFICATIONS

- Height: Min. 18.5" (± .5") 46.9cm (±1.27cm) at Seat Max. 37" at seat (93.9cm)
- Width: 27" (68.6cm) Length: 57" - 76" with foot section extended (144.7cm - 193cm)
- Leg rest: 21" W x 13" D (53.3cm W x 33cm D)
- Front Storage Drawer: 19" W x 12" D (48.26cm W x 30.48cm D)
- Rear Storage Drawer: 12.5" W x 24.5" D (31.75cm W x 62.23cm D)

**Brochure and images are a general overview of the product. Actual product options may vary.*