

# UMF 4070 ProGlide300 Power Exam Table



Stock Photo - Options May Vary

## FACTORY NEW

The UMF 4070 ProGlide300 Power Exam Table is designed for performance, patient comfort, and clinical efficiency. It features power Hi-Lo movement with an air-spring-assisted manual backrest and a one-function foot control for height adjustment.

The OneTouch WheelBase® central-locking caster system ensures safe, easy mobility. A low 18.5" ADA-compliant height allows better access, while the hand control offers two programmable positions and return-to-home.

The UMF 4070 ProGlide300 includes seamless cushions, dual storage drawers, four-position stirrups, and 7-inch pelvic tilt. It also features a stainless steel drain pan, drawer warmer, and ultra-comfort leg cushion with sturdy extension.

## FEATURES

- Power Hi-Lo with Air Spring Assisted Manual Backrest
- One-function Foot Control - Power Hi-Lo
- 650lb (294.84kg) Patient Weight Capacity
- OneTouch WheelBase® Integrated Central-Locking Caster System
- Low Access ADA Height 18.5" ( ±0.5") 46.9cm (±1.27cm)

## SPECIFICATIONS

- Height: Min. 18.5" ( ± .5") 46.9cm (±1.27cm) at Seat Max. 37" at seat (93.9cm)
- Width: 27" (68.6cm)
- Length: 57" - 76" with foot section extended (144.7cm - 193cm)
- Leg rest: 21" W x 13" D (53.3cm W x 33cm D)
- Front Storage Drawer: 19" W x 12" D (48.26cm W x 30.48cm D)

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)