

# Oakworks 3000 Series Procedure Chair



Stock Photo - Options May Vary

## FACTORY NEW

The Oakworks 3000 Series Procedure Chair is a fully powered, ADA-compliant solution designed for a wide range of clinical procedures. Its smooth, electric adjustments for height, backrest, leg rest, and Trendelenburg/reverse Trendelenburg enhance both patient comfort and provider efficiency. The hand control features programmable memory positions, allowing quick transitions between exam setups.

With an ADA-compliant low height of 19", patient transfers are safer and easier. The 2.5" Comfort Foam pad is removable for cleaning or replacement, while water-resistant fabric withstands common disinfectants. An open base design improves ergonomic access, and optional wide headrest extension adds further versatility.

Customize your table with 16 fabric color options. Includes an 18" paper roll holder.

## FEATURES

- Powered height, Trendelenburg/Reverse, Leg Rest, & Backrest adjustments
- Hand control has programmable memory positions to save time making chair adjustments
- ADA compliant 19" low height
- Highly water resistant fabric can be wiped down with disinfectants as needed
- 2.5" Comfort Foam removeable for cleaning/replacement

## SPECIFICATIONS

- Powered Height, Backrest, Trendelenburg (tilt), Leg Rest
- 550 lb bariatric, twin tower design
- Open base design for better access
- 2.5" Comfort Foam™ Padding
- 3 Section Top - 26in. (66cm) wide x 61 in. (155cm) long with Head Rest [76in. (193cm) long with Flat Head Rest and Platform extended

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)