

Oakworks 3100 Series Procedure Chair



Stock Photo - Options May Vary

FACTORY NEW

The OakWorks 3100 Series Procedure Chair is a premium, USA-made solution designed for imaging, pain management, and a wide range of medical procedures. With a rock-solid, fully powered four-motion system—height, backrest, Trendelenburg/reverse, and leg rest—this chair offers smooth, quiet positioning from 19" to 33" height range. Both hand and foot controls feature programmable memory positions for faster workflow.

The 3-section top provides ample surface space, extending to 76" with the flat headrest. The chair supports up to 550 lbs and includes retractable stirrups for patient comfort. It features an 18" paper roll holder for convenience.

Furthermore, A supine headrest system enhances patient positioning during exams.

FEATURES

- Motion 1 - Electric Height Range: 19in. - 33in. (48cm-84cm)
- Motion 2 - Electric Trendelenburg/Reverse Tilt: +15°/-21°
- Motion 3 - Electric Backrest: 0-75°
- Motion 4 - Electric Leg Rest: 0-90°
- Hand Control - Adjust Height, Trendelenburg/Reverse, Leg Rest, & Backrest with programmable memory positions

SPECIFICATIONS

- 3 Section Top - 26in. (66cm) wide x 61 in. (155cm) long with Head Rest [76in. (193cm) long with Flat Head Rest and Platform extended]
- Backrest - 25in. (63.5cm) L
- Seat - 24in. (61cm) L
- Leg rest - 11in. (28cm) L
- Motion 1 - Electric Height Range: 19" - 33" (48cm-84cm) low height and 22" - 36" (56-91cm) high height

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)