

Oakworks 3050 Series Procedure Chair



Stock Photo - Options May Vary

FACTORY NEW

OakWorks 3050 Series Procedure Chair. The Oakworks 3050 is a high-quality, U.S.-made chair built for comfort and easy use in medical and beauty settings.

Smooth electric controls and ADA-compliant low height make positioning easy and improve patient access on and off the chair. Adjustable armrests add flexibility, making it great for Med Spas, injectables, microneedling, and many other treatments.

This chair supports a robust 550-pound bariatric weight capacity, making it suitable for medical weight loss and wellness practices. Its sturdy build ensures it can handle diverse patient needs while maintaining safety and comfort.

FEATURES

- Motion 1 - Electric Height Range: 19in. - 33in. (48cm-84cm)
- Motion 2 - Electric Trendelenburg/Reverse Tilt: +15°/-21°
- Motion 3 - Electric Backrest: 0-75°
- Motion 4 - Electric Leg Rest: 0-90°
- Hand Control - Adjust Height, Trendelenburg/Reverse, Leg Rest, & Backrest with programmable memory positions

SPECIFICATIONS

- 3 Section Top - 26in. (66cm) wide x 64 in. (162cm) long with Head Rest [76in. (193cm) long with Flat Head Rest and Platform extended]
- Backrest - 25in. (63.5cm) L
- Seat - 24in. (61cm) L
- Leg rest - 11in. (28cm) L
- Electric Height Range: 19" - 33" (48cm-84cm) low height or 22" - 36" (56-91cm) high height

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)