

Hill-Rom Advanta 2 (P1190A) Hospital Bed - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Hill-Rom Advanta 2 (P1190A) Hospital Bed is designed to enhance comfort and support for both patients and caregivers. It features pressure redistribution technology and built-in guards to reduce patient sliding, helping prevent injuries and support patient safety. An integrated exit alarm aids in fall prevention efforts.

The bed includes nurse call integration that provides real-time data and alerts, enabling care teams to respond faster and more effectively. Positioning controls allow caregivers to place patients in safe, comfortable positions with minimal effort, contributing to faster recovery and improved patient experience.

The Advanta 2 is trusted in healthcare facilities worldwide as a reliable choice for providers seeking to enhance patient safety, comfort, and caregiver workflow. The bed includes a mattress.

FEATURES

- Four-corner brake and steer for easy mobility.
- One-button Dining Chair feature places patients at a low height.
- Built-in bed extender lets a single caregiver add length to the bed.
- Integrated Line Manager keeps IV lines untangled and out of the way.
- Care Grip patient egress assist handle.
- In-rail storage compartment.
- Siderail patient positioning controls are freely accessible to the patient.

SPECIFICATIONS

- Low position floor to deck: 15.75" (40.01 cm)
- High position floor to deck: 32.5" (82.55cm)
- Width: 36" (91.44 cm)
- Normal length: 84" (213.4 cm)
- Extended length: 88" (223.5 cm)
- Max head of bed elevation: 65°
- Max knee elevation: 36°
- Weight Limit (Safe working load): 500 lbs

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)