

Hill-Rom VersaCare P3200 Hospital Bed - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Hill-Rom VersaCare P3200 Hospital Bed delivers a high-performance solution for mid- to high-acuity patients. Its ergonomic design and advanced motorized functions improve patient positioning, reduce caregiver strain, and support better clinical outcomes.

Ideal for intensive care and post-surgical recovery, the VersaCare bed enhances both patient comfort and staff efficiency. Refurbished models provide healthcare facilities with a reliable, cost-effective alternative to new beds without sacrificing essential features or quality.

This Hill-Rom bed includes important safety features like bed exit alarms, weight sensors, and adjustable side rails. Moreover, these features help lower the risk of falls and make it easier for staff to care for patients safely.

FEATURES

- One-Button Boost feature helps patients move toward the head of the bed
- The Point-of-Care control panel uses easy to read digital displays and simple icons
- 3-Level Bed Exit Alarm with advanced sensitivity detects whether the patient is lying in bed, sitting on the edge of the bed, or is out of bed
- FlexAFoot platform extender customizes the bed length up to 12" and helps keep patients from sliding down the bed
- LowChair Function allows the patient to be placed in an upright sitting position while keeping the bed low

SPECIFICATIONS

- Expertly reconditioned
- Includes new foam mattress
- Weight Capacity: 500 lb.
- Sleep Surface Dimensions: 84" x 36"
- Low Position Height: 18"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)