

Hill-Rom Advanta (P1600) Hospital Bed - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Hill-Rom Advanta (P1600) Hospital Bed is fully electric, adjustable, and built for reliability and ease of use. It offers the ideal balance of function, comfort, and control in hospitals, training facilities, or home care settings.

The Advanta™ P1600 features a retractable frame and 4 motors for easy positioning, including Trendelenburg and chair-in/out. This flexibility enhances patient comfort and simplifies repositioning.

Furthermore, the bed prioritizes patient safety with features such as safer entry and exit, an optional Patient Positioning Monitoring (PPM) bed exit system, and fall-prevention alerts—helping caregivers reduce risks, monitor patient movement, and respond quickly to changes in patient status.

FEATURES

- Fully electric, adjustable hospital bed
- 4-motor, multi-function design for smooth operation
- Retractable bed frame for comfort and space efficiency
- Chair-in and chair-out positioning for easy patient access
- Trendelenburg and reverse Trendelenburg positioning

SPECIFICATIONS

- Low position (minimum height): 54" (137 cm)
- High position (maximum height): 74" (188 cm)
- Maximum width: 36" (91 cm)
- Total length: 91" (231 cm)
- Retracted length: 78.5" (199.4 cm)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)