

# Harloff Medium Narcotics Cabinet | NC12C16-ST2



*Stock Photo - Options May Vary*

## **FACTORY NEW**

The Harloff Medium Narcotics Cabinet NC12C16-ST2 provides secure storage for controlled substances in medical and pharmaceutical settings. Constructed from durable 20-gauge steel, it features a double-lock system with two high-security, pick-resistant tubular locks that require dual access to reduce unauthorized entry and help facilities maintain compliance.

The cabinet offers flexible interior organization with one adjustable or removable shelf and five mounting positions to accommodate various medication sizes and storage needs. Measuring 12" H x 16" W x 8" D and weighing 15 pounds, it fits efficiently into compact spaces without requiring extensive room.

Installation is straightforward with pre-drilled holes and included mounting hardware. The combination of reliable construction, dual-lock security, and convenient design makes this cabinet suitable for hospitals, clinics, and veterinary offices seeking to strengthen controlled substance protection.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)