

Harloff ADA-Compliant In-Room Medication Storage Cabinet | WV2795-CM



Stock Photo - Options May Vary

FACTORY NEW

The Harloff ADA-Compliant Medication Storage Cabinet WV2795-CM is designed for long-term care settings, combining secure storage with a slim profile and wood vinyl finish suitable for resident rooms. Constructed from heat-sealed powder-coated steel with an aluminum door, it provides strength and durability with an attractive wood-like appearance.

The cabinet holds up to 40 punch cards—20 per shelf—or a mix of punch cards and liquid medications using the included shelf and divided bin. A drop-down shelf provides added convenience during medication administration. Stainless steel pinned door hinges and a one-hand slam lock enhance both security and ease of use.

Each unit comes standard with a key lock and two keys. The cabinet can be upgraded to an electronic pushbutton lock for multi-user access. Backed by a 12-year manufacturer's warranty, the WV2795-CM offers reliable and accessible medication storage for care facilities.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)