

# Harloff Tall Narcotics Cabinet with Electronic Lock | NC24C16-SE1



Stock Photo - Options May Vary

## FACTORY NEW

The Harloff Tall Narcotics Cabinet NC24C16-SE1 is a heavy-duty wall-mounted cabinet constructed from 20-gauge steel, designed for secure storage of controlled substances in medical facilities. It features an electronic locking outer door with capacity for up to 20 user codes and field-reprogrammable lock functionality.

The cabinet measures 24" H x 16" W x 8" D and includes two adjustable or removable shelves that can be positioned at multiple heights to accommodate various medication bottles and containers. Predrilled holes and mounting hardware are provided for straightforward installation on walls or solid surfaces.

The NC24C16-SE1 meets DEA requirements for storing Schedule III, IV, and V controlled substances when placed within a locked office or secure enclosure. It is not suitable for Schedule I and II drugs without additional security measures. Consult current DEA guidelines to ensure compliance with applicable regulations.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)