

Bair Hugger 675 Patient Warming System – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Bair Hugger 675 Patient Warming System – Refurbished. The 3M Bair Hugger Warming Unit Model 675 provides effective forced-air warming with four temperature settings: Ambient, Low (89.6°F), Medium (100.4°F), and High (109.4°F).

It features an intuitive interface, ergonomic hose end design for easy blanket insertion, and alarms for over-temperature and under-temperature conditions, ensuring patient safety in all clinical settings. The Bair Hugger Model 675 is compatible with all 3M Bair Hugger Warming Blankets and Gowns, across all sizes.

In addition, it includes a detachable hose and MERV 14 filter with a change indicator. Additionally, the filter lasts up to 12 months or 500 hours.

FEATURES

- 4 temperature settings: 43°C, 38°C, 32°C, and ambient
- Small size frees up more usable work space
- Light weight makes transport and set-up a breeze
- Freestanding or easily attached to an IV pole, bedrail
- Modern look and compact design

SPECIFICATIONS

- Height: 14.5 in (36 cm)
- Width: 10.5 in (27 cm)
- Depth: 11 in (28 cm)
- Weight: 10 lbs (4.5 kg)
- High: 109.4°F (43°C)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)