

Bair Hugger 775 Patient Warming System – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Bair Hugger 775 Patient Warming System – Refurbished. The Bair Hugger 775 is a quiet, easy-to-use warming unit designed for effective patient temperature management. With adjustable airflow and precise temperature control, it delivers consistent warmth with minimal noise.

The unit features two airflow settings and hose-end temperature sensing, ensuring accurate and safe temperature delivery, tailored to each patient's needs. Equipped with three temperature sensors, the Bair Hugger 775 provides rapid response for quick heat transfer.

The adjustable airflow further enhances performance, offering optimal comfort and safety for patients of all sizes, from pediatric to adult. Additionally, it is compatible with various Bair Hugger blanket styles, ensuring flexibility for different clinical settings.

FEATURES

- Quiet, easy-to-use unit.
- Adjustable airflow, precise temperature delivery, and ease of use from a single, quiet warming unit.
- Safe, quick temperature response.
- 2 airflow settings and hose-end temperature sensing for precise temperature delivery.
- 3 temperature sensors and adjustable airflow provide safe, quick temperature response and greater heat transfer.

SPECIFICATIONS

- Height: 12" (30 cm)
- Width: 10" (25 cm)
- Depth: 13.5" (34 cm)
- Weight 15.5 lb (7 kg)
- High: $43^{\circ} \pm 1.5^{\circ}\text{C}$ ($109.4^{\circ} \pm 2.7^{\circ}\text{F}$)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)