

Pedigo Deluxe Combination Blanket & Fluid Warmer | P-210



Stock Photo - Options May Vary

FACTORY NEW

Pedigo P-210 Blanket/Fluid Warmer Cabinet. The Pedigo P-210 Blanket/Fluid Warmer Cabinet combines blanket and fluid warming in one efficient, deluxe unit. Its compact design fits easily in clinical spaces while supporting versatile warming needs.

The cabinet features a clear window door, allowing staff to quickly monitor inventory and maintain organized storage. This visibility improves workflow and reduces the need to open the door unnecessarily.

Equipped with a 4-digit LED display, the P-210 provides precise temperature control for both blankets and fluids. The control lock programming enhances safety, while the auto time programming simplifies scheduling and reduces human error.

FEATURES

- Ideal for countertop use
- The unit can be stacked
- Four digits LED display
- Control lock programming
- Auto time programming

SPECIFICATIONS

- Temperature: Up to 200° F
- Voltage: 125 VAC
- Compartment Capacity: 2.3 cubic feet
- Compartment Dimensions: 13¼ W x 21½"D x 13-13-7/8 H
- Overall Dimensions: 16-3/16 W x 25-7/8 D x 20-5/8 H

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)