

# Bair Hugger 750 Patient Warming System – Refurbished



Stock Photo – Options May Vary

## FULLY REFURBISHED • CLINICALLY CERTIFIED

Bair Hugger 750 Patient Warming System – Refurbished. The Bair Hugger 750 is designed for accurate, efficient temperature management. Hose-end temperature sensing ensures precise air temperature delivery, while three integrated temperature sensors guarantee rapid response.

The unit's increased airflow further enhances quick and effective warming, making it ideal for diverse clinical needs. This advanced system tracks key therapy metrics, including duration, over-temperature events, and calibration status.

Additionally, the built-in hour meter records usage, and fault codes are easily accessible via the front panel. There's no need to open the unit to access this essential information, streamlining operation and maintenance.

## FEATURES

- Hose-end temperature sensing ensures accurate air temperature delivery
- Three temperature sensors and increased airflow for quicker response
- Temperature settings track and report forced-air warming therapy duration
- Over-temperature, calibration, built-in hour meter and fault code reporting via the front panel – no need to open the unit
- Bair Hugger temperature management units are the ONLY forced-air warming systems that can be used with an integrated, cost-effective method of fluid warming – The Bair Hugger 241® blood/fluid warming set

## SPECIFICATIONS

- Size: 12.5" H x 13.5" D x 10" W, 31.5 cm H x 34 D x 25 cm W
- Weight: 16.3 lb (7.4 kg)
- 55 dBA
- Flexible, compatible with Bair Hugger brand 241 blood/fluid warming system.
- High efficiency 0.2 µm Filter

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)