

Bair Hugger 505 Patient Warming System – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Bair Hugger 505 Patient Warming System – Refurbished. The 3M Bair Hugger 505 sets the standard in forced-air warming, offering safe, quiet, and effective temperature control. Furthermore, its compact, lightweight design makes it perfect for busy healthcare environments, ensuring easy transport and quick setup.

Adaptable to any space, it works both freestanding or mounted on various supports, delivering reliable warmth from operating rooms to recovery wards. This versatile unit offers more than just patient warmth. It's the only forced-air system offering integrated fluid warming, providing a comprehensive, cost-effective temperature management solution.

Additionally, the Bair Hugger 505 operates quietly, providing optimal warmth without disrupting sensitive environments like operating rooms.

FEATURES

- Small size frees up more usable workspace
- Light weight makes transport and set-up a breeze
- Freestanding or easily attached to an IV pole, bedrail or optional rolling cart
- Built-in hour meter makes it easy to monitor usage for preventative maintenance
- Unique snap-fit hose swivels at three points for easy blanket attachment and positioning

SPECIFICATIONS

- 110-120 VAC, 60 Hz, 9.5 Amperes
- 220-240 VAC, 50 Hz, 4.5 Amperes
- 100 VAC, 50/60 Hz, 10 Amperes
- High-efficiency 0.2 µm air filter
- Meets UL 60601-1 requirements

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)