

Steris 5085 SRT Standard Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Steris 5085 SRT Standard Series Replacement Pad Set is a 3-piece cushion set designed for low-risk patients and short-duration surgical procedures. It provides reliable support for routine surgical use while offering a budget-conscious solution without compromising safety or essential comfort.

The pad set features 2-inch thick medium-firm high-density foam that delivers consistent support and stability. The electrically conductive vinyl cover helps reduce static and resists fluids and disinfectants, making it suitable for daily use in general surgical environments.

The lightweight construction allows for quick and easy pad replacement when needed. This standard series is ideal for hospitals where advanced pressure relief is not required for routine procedures.

FEATURES

- For Low-Risk Patients
- For Short-Duration Procedures
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex-Free
- Hook & Loop kit Included

SPECIFICATIONS

- Construction: Medium Firm High-Density Exterior
- Cover: Electrically Conductive Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)