

Surgical Armboard Replacement Pads



Stock Photo - Options May Vary

FACTORY NEW

Venture Medical offers replacement surgical armboard pads designed for operating table systems. These pads are available in multiple thickness options ranging from 1 inch to 4 inches to accommodate various table configurations and clinical requirements. Standard models measure 6 inches wide and come in 24-inch or 26-inch lengths.

All pads feature flame-retardant, anti-microbial, and fluid-resistant properties to meet surgical environment standards. The electrolyte cover is easy to clean and maintain for infection control. A specialized bariatric armboard pad (10 inches wide) is available for specific table models requiring extended width capacity.

Custom pad configurations are available upon request, with options for different foam styles, fabric thicknesses, and attachment devices. Standard models include loop Velcro down the center bottom unless otherwise specified at time of order.

FEATURES

- 1", 2", 3" or 4" thick to match table pad heights
- Easy to clean electrolyte cover
- Available 6" wide in 24" or 26" long
- Bariatric armboard pad is for item # MCM412

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)