

# MCM Round Hamper Stands



Stock Photo - Options May Vary

## FACTORY NEW

The MCM Round Hamper Stands deliver strength, hygiene, and reliability for demanding healthcare environments. Fully welded 16-gauge stainless steel construction enhances durability and stability, while the stainless steel surface is easy to clean and supports infection control. Smooth finishes reduce debris and allow faster sanitation.

The classic round design holds more while saving floor space, allowing staff to place the hamper easily in patient rooms, hallways, or utility areas. The versatile shape fits well into different layouts and daily workflows.

Mobility features improve efficiency during transport and collection. Three-inch swivel casters allow smooth movement across various flooring types, with full platform support keeping loads stable during transport. These stainless steel round hamper stands combine durability, hygiene, and ease of use for professional healthcare and laboratory operations.

## FEATURES

- Fully welded 16-gauge stainless steel construction
- Rust- and bacteria-resistant stainless steel surface
- Classic round design for maximum capacity
- Easy-to-clean and sanitize surfaces
- 3-inch swivel ball-bearing casters for smooth mobility
- Full platform support to prevent bag sagging
- Standard 35-inch height for ergonomic use
- Optional SayNo2Rust casters for enhanced corrosion resistance

## SPECIFICATIONS

- 18" or 25" Height
- Construction: Stainless Steel, Fully Welded 16-gauge
- Caster Diameter: 3"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)