

# MCM Triangle Hamper Stands



Stock Photo - Options May Vary

## FACTORY NEW

The MCM Triangle Hamper Stands provide practical, durable solutions for medical facilities, laboratories, and healthcare environments. Constructed from fully welded 16-gauge stainless steel, they withstand rigorous daily use while remaining easy to maintain. The stainless steel resists rust and bacteria, ensuring hygienic performance.

The unique triangle shape saves space, fitting neatly into corners or against walls, ideal for areas with limited floor space. Three-inch swivel wheels enable easy movement between rooms and down hallways. A solid platform keeps bags from sagging, making removal clean and convenient while preventing spills.

Available in two sizes with 18-inch and 25-inch height options, these hamper stands adapt to different facility needs. The sturdy design minimizes replacements, offering long-term cost-effective value. Stainless steel construction simplifies cleaning and sanitizing, supporting a safe, germ-free environment.

## FEATURES

- Fully welded 16-gauge stainless steel construction
- Unique triangle shape saves corner and wall space
- 3-inch swivel ball-bearing casters for smooth mobility
- Full platform prevents bag sagging and spills
- Rust and bacteria resistant, easy to clean
- Optional SayNo2Rust casters for corrosion-free use
- Durable, long-lasting design for cost-effective healthcare use

## SPECIFICATIONS

- Material: 16-gauge stainless steel
- Shape: Triangle, space-saving
- Height: 18 inches or 25 inches
- Casters: 3-inch swivel (optional SayNo2Rust)
- Rust and bacteria resistant

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)