

# MCM Lidded Triangle Hampers



Stock Photo - Options May Vary

## FACTORY NEW

The MCM Lidded Triangle Hampers are designed for hygiene, durability, and efficiency in medical, laboratory, and healthcare environments. The secure lid keeps waste and linen contained, supporting strict infection control standards. The unique triangle shape fits neatly into corners or against walls, saving valuable floor space while maintaining full capacity.

Staff can maneuver these hampers effortlessly thanks to 3-inch swivel ball-bearing casters, which provide smooth, stable movement in busy corridors and patient areas. A full platform design supports bags fully, preventing sagging and reducing spills, which enhances hygiene and simplifies bag removal.

These triangle hampers combine functionality with easy sanitation, resisting bacteria buildup and supporting efficient operations. Available in two sizes with a standard height, the fully welded stainless steel construction maximizes space, hygiene, and durability in any healthcare facility.

## FEATURES

- Fully welded 16-gauge stainless steel construction
- Unique triangle shape saves corner and wall space
- 3-inch swivel ball-bearing casters for smooth mobility
- Full platform prevents bag sagging and spills
- Rust and bacteria resistant, easy to clean
- Optional SayNo2Rust casters for corrosion-free use
- Durable, long-lasting design for cost-effective healthcare use

## SPECIFICATIONS

- Material: 16-gauge stainless steel
- Shape: Triangle, space-saving
- Height: 18 inches or 25 inches
- Sizes: Two options
- Casters: 3-inch swivel (optional SayNo2Rust)
- Rust and bacteria resistant

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)