

Precision Medical Air Chrome Flowmeter



Stock Photo - Options May Vary

FACTORY NEW

The Precision Medical Air Chrome Flowmeter has pressure-compensated performance for demanding clinical environments. This precision flowmeter meets strict global durability and accuracy standards. First, the chrome finish supports long-term use in high-traffic medical settings.

Then, the advanced internal design maintains consistent flow across low- and high-range applications. Additionally, the polycarbonate flow tube provides clear visibility and dependable readings. Clinicians gain precise control during sensitive procedures.

The needle valve design supports accurate flow adjustments with confidence. As a result, users achieve stable output for neonatal, pediatric, and adult care. The stork icon clearly identifies the low-flow neonatal model at a glance.

FEATURES

- Highly accurate needle valves ensure precise flow-settings
- Supply pressure of 50 psi
- Air and oxygen models available
- Yellow knob marked with "AIR" for easy identification
- Durable polycarbonate housing

SPECIFICATIONS

- Model: 1MFA1XXX
- Flow Range: 0-15 l/min
- Graduations: 0.5 l/min (0-5), 1.0 l/min (5-15)
- Accuracy
- 0-5 l/min \pm 0.25 l/min
- 6-15 l/min \pm 0.50 l/min
- Model: 1MFA11XX
- Flow Range: 2-26 l/min

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)