

# Vapofil Vaporizer Filler Adapters

## FACTORY NEW

Vapofil Vaporizer Filler Adapters simplify anesthetic filling for key-filled vaporizers by replacing traditional filler keys with a faster, cleaner, and more controlled filling solution. The dual-tube design features one tube for air escape and a second tube for smooth anesthetic agent delivery, actively preventing vapor-lock and ensuring consistent flow while minimizing waste during refills.



Stock Photo - Options May Vary

The compact construction supports convenient handling in busy medical settings while maintaining reliable, accurate agent transfer. Healthcare professionals trust Vapofil for its secure design and ease of use, reducing downtime and keeping procedures moving smoothly.

Vapofil is available in two color-coded options: Ultane/Sevoflurane (yellow) and Forane/Isoflurane (purple) for quick identification and improved safety organization. The adapter meets CE compliance standards and includes California Proposition 65 compliant labeling.

## FEATURES

- Designed for key-filled anesthetic vaporizers
- Dual inner tube system prevents vapor-lock
- Allows smooth anesthetic agent flow
- Enables air to escape during filling
- Reduces spills and agent waste
- Improves filling speed and efficiency
- Supports cleaner, controlled transfers
- Color-coded options for quick identification

## SPECIFICATIONS

- Application: Anesthesia vaporizer filling
- Compatibility: Key-filled vaporizers
- 2 Available Types: Ultane/Sevoflurane (Yellow), and Forane/Isoflurane (Purple)
- Compliance: CE Certified
- MFG P/N: 8907-S (Sevo/Ultane)
- MFG P/N: 8907-F (Iso/Forane)
- Regulatory Warning: California Prop 65 compliant labeling included

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)