

Allen Yellofin (Yellowfin) Lift Assist Stirrups – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Allen Yellofin Lift Assist Stirrups enable surgical personnel to lift and position a patient's legs, then safely and securely hold them in position during surgical procedures. The stirrups feature a lateral safety fin design that eliminates the possibility of perineal nerve injury, unlike candy cane stirrups that can cause nerve damage from uncontrolled abduction.

The patented Lift-Assist technology neutralizes leg weight, allowing staff to move the leg to any desired position easily. This is particularly helpful with heavy or obese patients, preventing staff back injuries while enhancing patient comfort and safety. The stirrups attach to any OR table and come with new pads and table clamps.

FEATURES

- Lateral safety fin eliminated the possibility of perineal nerve injury
- Patented Lift-Assist™ technology Lift-Assist technology neutralizes the leg weight
- Easy to position the leg to any desired position, even heavy/obese patients
- Aids in preventing staff back injuries
- Enhances patient comfort and safety
- Attaches to any OR table
- Includes: New Pads and Table Clamps

SPECIFICATIONS

- Lithotomy Range: Any position between +84° to -33° lithotomy
- Abduction Range: Set the stirrup in any position between +25° to -9° abduction
- Patient Weight Limit: 350 lb. patient (159 kg)
- Includes: New Pads and Table Clamps

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)