

# Drager Fill Vaporizer Filler Adapters



Stock Photo - Options May Vary

## FACTORY NEW

Drager Fill Vaporizer Filler Adapters are designed to improve accuracy and safety when filling anesthesia vaporizers. These adapters support precise control during agent filling procedures and help reduce spills, maintaining proper vaporizer function. Made from durable materials with consistent fit, they perform reliably in busy clinical settings.

Two clearly identified adapter types are available to match your anesthetic agent. The Ultane/Sevoflurane adapter (Yellow, SKU M36120) supports sevoflurane applications, while the Forane/Isoflurane adapter (Purple, SKU M36110) supports isoflurane use. Color coding enables quick selection of the correct option for your setup.

Each adapter installs easily and creates a tight, secure connection during filling. The design supports smooth workflow and accurate agent transfer, helping staff work efficiently while ensuring safe and controlled anesthesia delivery in surgical and clinical settings.

## FEATURES

- Designed for use with Drager Fill vaporizers for precise anesthesia agent filling
- Two adapter types clearly color-coded for easy identification
- Helps reduce anesthetic spillage and improves filling accuracy
- Durable construction for reliable performance in clinical settings
- Quick and easy to install during agent refill procedures
- Supports safe and efficient anesthesia control for clinicians

## SPECIFICATIONS

- Compatibility: Works with Drager Fill vaporizer systems
- Adapter Types: • Ultane/Sevoflurane (Yellow) - SKU M36120 • Forane/Isoflurane (Purple) - SKU M36110
- Material: High-quality medical-grade components
- Color Coding: Yellow for sevoflurane agents, purple for isoflurane agents
- Application: Anesthesia vaporizer agent filling
- Use Case: Clinical and surgical anesthesia delivery environments

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)