

Luvis S200 LED Minor Surgical Light Series



Stock Photo - Options May Vary

FACTORY NEW

The Luvis S200 LED Minor Surgical Light Series is perfect for minor procedures and examinations in small spaces. The S200 provides a shadowless effect that offers less shadow, even under a 50% screened condition.

The fanless design generates no heat or noise. Are you looking for a specific LED Surgical Light or additional types of surgical lights? Click any of the following links to view our complete catalog of LED Surgical Lights, Surgical/OR Lights, Minor Surgical/Procedural Lights, or Exam/Diagnostic Lights.

Have Questions? Contact us today for questions or available options to fit your specific needs!

FEATURES

- Excellent Illumination
- Ergonomic & Friendly Design
- Detachable & Autoclavable Main Handle
- No UV Emission & Cool Light
- Excellent Depth of Illumination

SPECIFICATIONS

- LEDs: 18
- CRI: 90
- Color Temperature: 3,800K /4,300K/4,800K
- Head Size: 13.78" x 13.39"
- Focus Size: 7" (Circle)/11.8" (Circle)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)