

# Burton Wave LED Magnifier | WML35W45EC120



*Stock Photo - Options May Vary*

## FACTORY NEW

The Burton Wave LED Magnifier delivers clear, shadow-free visibility for detailed examinations in medical, dermatology, and laboratory environments. It combines a 3.5-diopter optical glass lens with dual independently controlled LED modules for customizable lighting. Users can switch between symmetrical lighting for shadow-free viewing and asymmetrical lighting for enhanced 3-D depth perception.

The magnifier features a self-balancing arm and hands-free adjustment capability, allowing users to position the light easily for optimal performance. Its ergonomic design integrates seamlessly into clinical settings while supporting long-term reliability with energy-saving features including automatic shut-off and step dimming.

The articulating K-arm design provides complete freedom of movement without drift, and a standard table edge clamp mount is included for secure positioning. The rectangular optical-quality glass lens ensures sharp zoom with consistent performance and superior visibility.

## FEATURES

- Shadow-free magnification using symmetrical lighting or 3-D magnification using asymmetrical magnification
- Self-balancing shade and hands-free neck assembly make precise adjustments easy
- Articulating K-arm standard
- Rectangular white crown optical-quality glass lens
- Equipped with two 6W individually controllable LED modules
- Automatic 9/4 hour auto shut-off feature and (0-50-100%) step dimming to support energy savings
- No drift spring-balanced arm allows complete freedom of movement
- Standard table edge clamp mount included

## SPECIFICATIONS

- Articulating Arm Length: 45"
- Glass Lens: 6.75" x 4.5" rectangular white crown optical-quality
- Optics: 3.5 (1.88X) Diopter
- Focal Length: 11"
- Illumination of 4600 lux
- Color Temperature of 4000 K
- CRI (Color Rendering Index) of 80
- Average LED Life: 50,000 Hours

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)