


# SchureMed Berchtold Operon D850, D750 Surgical Table Replacement Pads



*Stock Photo - Options May Vary*

## FACTORY NEW

SchureMed replacement pads are designed for Berchtold Operon D850/D750 surgical tables and feature a 20" leg extension. These pads utilize multiple layers of engineered foam to maximize pressure-reducing properties and conform to the patient's body while preventing pressure on bony prominences.

The multi-layer foam construction combines different foam types to achieve optimal pressure management. No single foam material can effectively conform to the patient, provide adequate support, and prevent bottoming out simultaneously, making the layered approach essential for surgical patient comfort and safety.

SchureMed brings over 20 years of experience in pressure management solutions for operating room tables, stretchers, arm boards, and other surgical padding applications. Seven pad options are available to meet specific pressure management objectives for your facility.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)