

SchureMed CamLoc Clamp | 800-0271



Stock Photo - Options May Vary

FACTORY NEW

The SchureMed CamLoc Clamp (800-0271) is designed for use in the sterile field and attaches to operating table rails at any location, accommodating multiple drapes. It accepts accessories with round posts up to 5/8" diameter and flat mounting positions to 1" x 3/8" diameter.

The clamp supports a patient weight capacity of 500 lbs and is compatible with most OR table designs and brands. Installation is straightforward and can be performed by a single person without specialized tools or maintenance.

The device is latex-free, maintenance-free, and backed by a two-year warranty covering manufacturing defects. Under normal use conditions, the clamp has a five-year service life.

FEATURES

- Product not made with (NRL) Natural Rubber Latex
- The device supports a 500 lb. proportional patient load (6'4" tall patient per 99% human body model)
- The product warranty covers the product from manufacturing defects for two years
- Product is maintenance-free; check product condition before subsequent use
- The life of the device is five years under normal use
- Store device between -4 °F to +86°F

SPECIFICATIONS

- Length: 4 3/4" +/- 0.5"
- Width: 5" +/- 0.5"
- Depth: 2 1/4" +/- 0.5"
- Device Weight: 2.0 +/- 0.5 lbs.
- Attaches to the rail of the surgical table at any point on the rail
- Single-person installation
- Side rail dimensions: 0.374" X 1.122" (9.5 mm x 28.5 mm)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)