

SchureMed Deluxe Rail Clamp | 800-0248



Stock Photo - Options May Vary

FACTORY NEW

The SchureMed Deluxe Rail Clamp (800-0248) is a surgical table accessory designed to securely mount equipment and positioning devices to operating room tables. It accommodates both 1" x 3/8" flat mounting posts and round mounting posts up to 5/8" diameter, allowing flexible positioning of accessories along the table rail.

This clamp is compatible with most OR table designs and brands, attaching at rail ends or notches. Accessories can be mounted vertically or horizontally to suit specific surgical requirements across gynecological, urological, laparoscopic, and orthopedic procedures.

The device is maintenance-free and supports a patient weight capacity of 500 lbs. Installation is straightforward and can be completed by a single person at any point along the surgical table rail.

FEATURES

- Designed for a wide range of purposes; affixing most operating room table accessories to most OR table designs and brands
- Attaches at OR table rail ends or notches
- Accessories can be mounted either vertically or horizontally
- Patient weight capacity is 500 lbs (227 kg)
- The clamp is maintenance-free
- The life of the device is five years under normal use
- Attaches to the rail of a surgical table at any point on the rail
- Single-person installation

SPECIFICATIONS

- Length: 3" +/- 0.5" Handle Collapsed
- Width: 2" +/- 0.5" Handle Extended
- Depth: 4" +/- 0.5" (with rail handle)
- Device Weight: 1.00 +/- 0.5 lbs
- Holds accessory equipment up to 1 x 3/8" flat mounting post and round mounting posts up to 5/8"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)