

SchureMed Simple Clamp | 800-0228



Stock Photo - Options May Vary

FACTORY NEW

The SchureMed Simple Clamp (800-0228) is a surgical table accessory clamp designed to securely lock and position surgical equipment on operating room tables. It features both vertically and horizontally mounting blade options, allowing flexible placement of accessories anywhere along the OR table rails or at rail ends and notches.

The clamp accommodates flat mounting blades measuring 1" x 3/8"D for vertical mounting and 1" x 1/4"D for horizontal mounting, making it compatible with most standard OR table designs and brands. Installation is straightforward and can be completed by a single person without specialized tools or maintenance.

The device is built for durability and reliability in the demanding surgical environment, with a five-year lifespan under normal use conditions and a two-year warranty against manufacturing defects.

FEATURES

- Simple Clamp securely locks surgical accessories with flat 1" x 3/8"D vertically mounting blades
- Mounts anywhere along the operating room table rails 1" x 1/4"D horizontally mounting blades
- For most operating room table accessories to most OR table designs and brands
- Attaches at operating table rail end or notches
- Surgical accessories can be mounted either vertically or horizontally
- Product is maintenance-free
- The life of the device is five years under normal use
- Single-person installation

SPECIFICATIONS

- Length: 3.00" +/- 0.5"
- Width: 1.75" +/- 0.5"
- Depth: 1.375" +/- 0.5"
- Device Weight: 1.00 +/- 0.5 lbs.
- Holds accessory equipment with 1" x 3/8" flat vertically mounting blades
- For Accessories with 1" x 1/4"D (2.5 cm x .63 cm) horizontally mounting blades
- Patient weight capacity: 720 lbs
- Clamp warranty covers the product from manufacturing defects for two years

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)