

SchureMed Split Leg Clamp | 800-0336



Stock Photo - Options May Vary

FACTORY NEW

The SchureMed Split Leg Clamp (800-0336) is a surgical table accessory designed to securely lock and position equipment along operating room table rails. It features mounting blades measuring 1" x 3/8"D (2.5 cm x .95 cm) for vertical mounting and 1" x 1/4"D (2.5 cm x .63 cm) for horizontal mounting, allowing flexible positioning of surgical accessories.

The clamp affixes most operating room table accessories and attaches at the operating table rail end or notches. SchureMed Clark Sockets work with various table designs to accommodate a wide range of surgical equipment and positioning needs.

Surgical accessories can be mounted either vertically or horizontally, providing adaptability for gynecological, urological, laparoscopic, and orthopedic surgical procedures. The clamp supports a patient weight capacity of 720 lbs (327 kg) and requires single-person installation.

FEATURES

- Securely locks surgical accessories with flat 1" x 3/8"D (2.5 cm x .95 cm) vertically mounting blades anywhere along the operating room table rails. 1" x 1/4"D (2.5 cm x .63 cm) horizontally mounting blades
- Clark Sockets by SchureMed are designed for a wide range of purposes
- Clamp affixes most operating room table accessories
- Attaches at operating table rail end or notches
- Surgical accessories can be mounted either vertically or horizontally

SPECIFICATIONS

- The patient's weight capacity is 720 lbs (327 kg)
- Single-person installation
- The product warranty covers clamp from manufacturing defects for 2 years
- The clamp is maintenance-free

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)