

SchureMed Lithotomy Stirrups | 800-0049



Stock Photo - Options May Vary

FACTORY NEW

The SchureMed Lithotomy Stirrups (Model 800-0049) provide secure patient positioning during lithotomy procedures. The extended lateral boot design fully supports the patient's leg, preventing superficial nerve damage and unwanted leg rotation during surgery. Extended soft boot liners are included for patient comfort.

These stirrups are suitable for both long and short procedures and accommodate patients up to 350 lb. The stirrups are sold in pairs and require two Schure Socket Clamps (sold separately) for attachment to OR table side rails.

Proper installation requires an OR table with side rails extending from the back to the foot section on each side.

FEATURES

- Use for long or short procedures
- Extended soft boot liners are included
- Weight capacity: 350 lb. patient weight
- Lithotomy Stirrups: Sold in pairs only
- Two required Schure Socket Clamps sold separately

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)