

SchureMed Pediatric E-Z Lift Stirrups | 800-0342-KID



Stock Photo - Options May Vary

FACTORY NEW

The SchureMed Pediatric E-Z Lift Stirrups (model 800-0342-KID) are designed for pediatric surgical procedures with a focus on ease of use and patient positioning. The stirrups feature exclusive axial rotation from the handle to ensure superior patient positioning and optimal support. The stirrup boots adjust easily with a single clamp, simplifying the process of positioning young patients, while included boot pads enhance comfort and provide a secure fit.

These stirrups offer improved stability and ease of use for surgical teams, helping reduce setup time while ensuring patient safety and comfort throughout the procedure. The floating boot automatically adjusts when repositioning the stirrup, and the design allows for easy control of the patient's leg when moving into operative position. Two SpringLoc Clamps are required for operation.

FEATURES

- Easy-to-use
- Easy to control patient leg when moving into operative position
- Exclusive axial rotation from the handle for superior patient positioning
- Boot adjusts with a single clamp which simplifies patient positioning
- Floating boot automatically adjusts when repositioning the pediatric boot stirrup

SPECIFICATIONS

- Set of 2
- Abduction +25°
- Adduction -9°
- Lithotomy range of +84° to -35° & externally rotates out 20°
- Weight Capacity: 250 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)