

SchureMed TKR Positioner for Total Knee Replacement | 800-0141 or 800-0171



Stock Photo - Options May Vary

FACTORY NEW

The SchureMed TKR Positioner is a ball-track system designed for total knee replacement procedures. The ball joint at the heel of the leg holder allows the operative leg to be positioned in extension, flexion, tilt, and rotation. Once locked, the system is completely rigid and easy to use.

This positioner eliminates the need for additional leg positioning devices, sandbags, and gel pads, streamlining the surgical setup. All components are autoclavable for sterilization and reuse.

The system requires a SchureMed CamLoc Clamp for installation on the surgical table.

FEATURES

- Track Locking Bracket securely locks the Base Plate to the table with two thumb screws
- Locks in any position with a simple turn of the wrist
- Smooth sliding motion for effortless leg positioning
- Well leg lies flat - eliminating the need for well-leg positioners
- All components are autoclavable

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)