

# SchureMed Bi-Lateral Leg Prepper | 800-0236



Stock Photo - Options May Vary

## FACTORY NEW

The SchureMed Bi-Lateral Leg Prepper model 800-0236 is designed to safely and effectively hold both legs during bilateral surgical preparations. The device features holders with a unique rotating and pivoting design that allows for unlimited configurations of the lower extremity to accommodate various surgical positioning needs.

The unit comes equipped with 2" Deluxe Pads for patient comfort and protection. The pads are removable, enabling straightforward cleaning and maintenance between procedures.

This leg prepper requires a Schure Socket XL for installation, which is sold separately.

## FEATURES

- Pads remove for easy cleaning
- Required Schure Socket XL

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)