

Clinton Powder Board Table | 7423

FACTORY NEW

Clinton Powder Board Table: 7423 (Previously 75-23K). Clinton Physical Therapy Equipment offers one of the best values in its class. It is tough for everyday physical therapy use and affordably priced for nearly every practice and facility.

Clinton offers a full line of treatment tables, staircases, mat platforms, parallel bars, and more. Their physical therapy equipment meets a wide range of needs. All of this comes with one source of convenience.

Clinton's Powder Board Table features an All Laminate Top in Maple or Gray, Oval Shape with Patient Cutout, Height Adjustable with Ultra-Heavy, Square, Telescoping Uprights.



Stock Photo - Options May Vary

FEATURES

- All laminate tops in Maple or Gray
- Tee molding edge bumper (Maple has tan tee molding, Gray has black tee molding)
- Oval shape with patient cutout
- White high-pressure laminate top with tee-molding edges
- two knobs adjust the height

SPECIFICATIONS

- Model Number: 7423
- Length: 63"
- Width: 40"
- Height Range: 26" to 34"
- Laminate: 1-1/8" thick

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)