

Clinton Personal Work Table w/ Small Cut-Out and Tilt-Top | 7433

FACTORY NEW

Clinton Personal Work Table w/ Small Cut-Out and Tilt-Top: 7433 (previously 76-33K). Clinton Physical Therapy Equipment offers one of the best values in its class.



Stock Photo - Options May Vary

It withstands the rigors of everyday physical therapy use and offers an affordable price for nearly every practice and facility. Clinton offers a full line of treatment tables, staircases, mat platforms, parallel bars, and more.

They design their equipment to meet a wide range of needs. All of this offers one source of convenience.

FEATURES

- Top tilts from 0° —50°
- All laminate top in Maple or Gray
- T-molding edge bumper (Maple have tan T-moulding, Gray has black T-moulding)
- Height adjustable legs with self-indexing, spring-loaded plungers
- Easy read numbered height indicator strip (no holes to count)

SPECIFICATIONS

- Model Number: 7433
- Length: 40"
- Width: 24"
- Height Range: 23" to 33"
- Top tilts from 0° —50°

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)