

# Clinton U-Shaped Table | 7425

## FACTORY NEW

Clinton U-Shaped Table: 7425 (previously 75-25K). Clinton Physical Therapy Equipment offers one of the best values in its class. It is built tough for everyday physical therapy use and is affordably priced for nearly every practice and facility.

Clinton offers a full line of treatment tables, staircases, mat platforms, parallel bars, and more. They provide a variety of physical therapy equipment to meet a wide range of needs.

All of this comes with one source of convenience. Clinton's U-Shaped Table features All laminated top in your choice of Maple or Gray, Height Adjustable legs with self-indexing, spring-loaded plungers.



*Stock Photo - Options May Vary*

## FEATURES

- All laminate top in Maple or Gray
- T-molding edge bumper (Maple has tan T-moulding, Gray has black T-moulding)
- Height adjustable legs with self-indexing, spring-loaded plungers
- Easy read numbered height indicator strip (no holes to count)
- The shape allows 4 wheelchair access

## SPECIFICATIONS

- Model Number: 7425
- Length: 72"
- Depth: 48" overall
- Height Range: 22" to 33"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)