

# Steris Harmony LED 585 Surgical Lights - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

Steris Harmony LED 585 Surgical Lights - Refurbished. The Steris Harmony 585 LED surgical light gives bright, clear lighting for almost any surgical setting. The 24" light head moves easily and uses long-lasting LEDs with a 30,000-hour life.

Each light comes fully refurbished, offering strong performance at a lower price. A 6-month warranty is included, making this a smart and affordable choice for hospitals, surgery centers, and medical clinics. In addition, the Steris 585 offers cutting-edge LED technology designed for modern surgical and clinical environments.

Furthermore, it delivers pure white, natural light with a 4400K color temperature and 14,800 ft candles (160,000 lux) intensity.

## FEATURES

- Cutting Edge LED Technology
- Pure White, Natural Light
- Excellent Depth of Illumination (26")
- Color temperature 4400k
- 360° Spring/Swivel Arm Support

## SPECIFICATIONS

- Intensity Pattern (large): 11"
- Intensity Pattern (small): 7"
- Clearance (rotation): 16' (approx)
- LED Life: Approx 30,000 hrs
- Intensity: Up to 14,800 Ft Candles

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)