

Luvis L400 Laminar Flow LED Surgical Light | Dual Ceiling Mounted



Stock Photo - Options May Vary

FACTORY NEW

The Luvis L400 Dual Ceiling Laminar Flow LED Surgical Light System is perfect for any procedure. The L400 features exceptional quality with an energy-efficient operation.

In addition, the L400 features an open design, which allows heat that may normally be “trapped” by the light head surface to escape. Additionally, the 96 LEDs in each light head produce up to 160,000 lux in a pattern that can be adjusted from 7.8” to 11.8”. The Luvis Surgical Light also features a control panel synchronized with a sub-handle, allowing the user to operate it easily during surgery.

Furthermore, smart Sensors automatically detect obstacles in the light field and increase intensity to compensate for reduced light.

FEATURES

- Excellent Illumination
- Excellent Dilution Effect
- Ergonomic & Friendly Design
- Detachable & Autoclavable Main Handle
- No UV Emission & Cool Light

SPECIFICATIONS

- LEDs: 96
- CRI: 97
- Color Temperature: 3,800k/4,300k/4,800k/5,400k
- Head Size: 28.75”
- Focus Size: 7” (Circle) - 11.8” (Circle)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)