

Luvis M300 LED Surgical Light | Dual Ceiling Mounted



Stock Photo - Options May Vary

FACTORY NEW

The Luvis M300 Dual Ceiling Mounted LED Surgical Light System is perfect for any surgical procedure. Additionally, the M300 features exceptional quality with an energy-efficient operation.

In addition, shop Luvis M300 for advanced dual LED surgical lighting. Improve clinical accuracy and bedside care with this high-intensity light.

Buy at Venture Medical. The Luvis Surgical Light also features: Head-mounted control panel Cool light with no UV emissions High Color Rendering Index of 95 Smart Sensors automatically detect obstacles in the light field and increase the intensity to compensate for a reduction in light.

FEATURES

- Excellent Illumination
- Excellent Dilution Effect & Ergonomic-Friendly Design
- Detachable & Autoclavable Main Handle
- No UV Emission & Cool Light
- Excellent Depth of Illumination

SPECIFICATIONS

- LEDs: 84
- CRI: 95
- Color Temperature: 3,800k/4,300k/4,800k
- Head Size: 23.6"
- Focus Size: 7.8" (Circle)/11.8" (Circle)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)