

Luvis M200 LED Surgical Light | Dual Ceiling Mounted



Stock Photo - Options May Vary

FACTORY NEW

The Luvis M200 Dual Ceiling Mounted LED Surgical Light System delivers exceptional illumination for a wide range of surgical procedures. The system features 54 LEDs in each lighthead that produce up to 120,000 lux with an adjustable illumination pattern ranging from 7.8" to 11.8" diameter. The 4,300k color temperature provides optimal visibility without UV emission or heat generation.

The M200 incorporates an ergonomic design with a detachable and autoclavable main handle for easy sterilization. A synchronized control panel integrated with the sub-handle allows surgeons convenient access to intensity and focus adjustments during procedures. The lightweight suspension system supports smooth positioning and stability.

The system is available in single, dual, or mobile configurations to accommodate various operating room layouts. It operates on standard 100-240V AC power with low energy consumption of 85W and includes a 2-year manufacturer's warranty.

FEATURES

- Detachable & Autoclavable Main Handle
- No UV Emission & Cool Light
- Intensity: 120,000 lux
- Color Rendering Index: 95
- Focusable Light Head

SPECIFICATIONS

- LEDs: 54
- Intensity: Max. 120,000lux
- Focus Size: 7.8" (Circle)/11.8" (Circle)
- Color Temperature: 3,800k/4,300k/4,800k
- LED Life Span: 50,000hr

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)