

# Medical Illumination System Two LED Duo w/ Monitor Arm | XLDS-S2MA



*Stock Photo - Options May Vary*

## FACTORY NEW

Medical Illumination System Two LED Duo w/ Monitor Arm: XLDS-S2MA. The energy-efficient System Two LED consumes half as much electricity as its halogen counterpart. In addition, long-life LEDs last 50,000 hours and reduce the need to purchase and store replacement bulbs.

They also eliminate procedure room downtime caused by unexpected bulb burnouts. Additionally, the XLDS-S2MA delivers 120,000 lux at a 4300K color temperature and reduces shadows for clearer visibility. Its MI arm system provides precise positioning and optimal surgical site access.

Furthermore, using System Two LED cool technology improves patient and staff comfort by reducing the heat produced by old-fashioned Halogen lights.

## FEATURES

- Green Technology
- Consumes half the electricity of its halogen counterpart
- Long-life 50,000-hour rated LEDs
- Output rated at 120,000 Lux with a Kelvin color temperature of 4300°
- Superior shadow reduction and designed

## SPECIFICATIONS

- Light Head Dimensions: 25.2" Diameter x 4.75" Deep
- Arm (Articulated): 4.0" Diameter (tapered) x 26.0" Long
- Ceiling Casting: 17.0" Diameter x 4.7" Deep
- Ceiling Rod: 2.5" Diameter x 9.0"-43.0" Long (depending on ceiling height)
- Duo Ceiling Mount: Approx. 104.5 lbs

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)